

Rabbit Anti-C3Ar1 Polyclonal Antibody For Wb, Elisa - Q16581

Item Number: CM0023058



Introduction

Rabbit polyclonal antibody against human C3AR1, a G protein-coupled receptor for anaphylatoxin C3a involved in chemotaxis and inflammation. Validated for WB and ELISA; cross-reacts with mouse. Useful for studying GPCR signaling.

[Learn More](#)

Parameter	Value
Product Name	C3AR1 Rabbit pAb
Remarks/Alias	AZ3B; C3AR; HNFAG09; C3AR1
Species	Human
Gene ID (Human)	719
Gene ID	719
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human C3AR1 (NP_004045.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q16581
Protein Weight	54kDa
Shipping	Ice bag