

Rabbit Anti-Flotillin-2 (Flot2) Polyclonal Antibody For Wb And Elisa - Q14254

Item Number: CM0021338



Introduction

Rabbit polyclonal antibody targeting human Flotillin-2 (FLOT2), validated for WB and ELISA. Cross-reacts with mouse. Suitable for studies on caveolae, cell adhesion, and epidermal biology.

[Learn More](#)

Parameter	Value
Product Name	FLOT2 Rabbit pAb
Remarks/Alias	ESA; ECS1; ESA1; ECS-1; M17S1; FLOT2
Species	Human
Gene ID (Human)	2319
Gene ID	2319
Immunogen	Recombinant protein corresponding to a sequence within amino acids 239-428 of human FLOT2 (NP_004466.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q14254
Protein Weight	47kDa
Shipping	Ice bag