


Rabbit Anti-Gca Polyclonal Antibody For Wb, If/Icc, Elisa - P28676

Item Number: CM0028155



Introduction

Rabbit anti-GCA polyclonal antibody targets calcium-binding grancalcin (~24 kDa). Applications: WB, IF/ICC, ELISA. Reactivity: human, mouse. Ideal for neutrophil adhesion and focal adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	GCA Rabbit pAb
Remarks/Alias	GCL; GCA
Species	Human
Gene ID	25801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-217 of human GCA (NP_036330.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28676
Protein Weight	24kDa
Shipping	Ice bag