

Rabbit Anti-Human IgG2 (Fc) Monoclonal Antibody For Wb, Elisa - P01859

Item Number: CM0010658



Introduction

Rabbit monoclonal antibody against human IgG2 (Fc) region, validated for WB and ELISA. Specific to human IgG2 Fc. Suitable for humoral immunity and antibody isotype research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit Anti-Human IgG2 (Fc) mAb
Species	Human
GeneID (Human)	3501
GeneID	3501
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-326 of human IgG2 (Fc) (P01859)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01859
Protein Weight	44kDa
Shipping	Ice bag