

Rabbit Anti-Mouse Cd117/C-Kit Monoclonal Antibody For Flow Cytometry - P05532

Item Number: CM0009637



Introduction

-conjugated rabbit anti-mouse CD117/c-Kit monoclonal antibody, validated for flow cytometry. Reacts specifically with mouse CD117, ideal for hematopoietic stem cell, mast cell, and melanogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD117/c-Kit mAb
Remarks/Alias	W; Bs; Fdc; Ssm; SC01; SC05; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3
Species	Mouse
Gene ID (Human)	3815
Gene ID	16590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-523 of mouse CD117/c-Kit (NP_001116205.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05532
Protein Weight	22kDa/108kDa/109kDa
Shipping	Ice bag