

Robo4 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wz75

Item Number: CM0024947

Introduction

ROBO4 Rabbit Polyclonal Antibody validated for WB and ELISA applications. Reacts with human ROBO4, immunogen corresponds to amino acids 848-1007. Expressed in endothelial cells, involved in angiogenesis and vascular patterning. Ideal for vascular biology and endothelial barrier studies.

[Learn More](#)

Parameter	Value
Product Name	ROBO4 Rabbit pAb
Remarks/Alias	MRB; AOVD3; ECSM4; ROBO4
Species	Human
GeneID (Human)	54538
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 848-1007 of human ROBO4 (NP_061928.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WZ75
Protein Weight	107kDa
Shipping	Ice bag