


Samm50 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Y512

Item Number: CM0026418



Introduction

Rabbit polyclonal antibody targeting human SAMM50 (mitochondrial cristae maintenance). Validated for WB, IHC-P, IP, ELISA. Cross-reacts with mouse, rat. Ideal for respiratory chain assembly research.

[Learn More](#)

Parameter	Value
Product Name	SAMM50 Rabbit pAb
Remarks/Alias	OMP85; SAM50; TOB55; TRG-3; CGI-51; YNL026W; SAMM50
Species	Human
GeneID (Human)	25813
GeneID	25813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SAMM50 (NP_056195.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y512
Protein Weight	52kDa
Shipping	Ice bag