


Sars-Cov-2 Nsp10 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P0Dtd1

Item Number: CM0023794



Introduction

Rabbit polyclonal antibody targeting SARS-CoV-2 NSP10. Suitable for Western blot, immunoprecipitation, and ELISA. Reacts with SARS-CoV-2. Ideal for studying viral replication, cap methylation, and host immune evasion.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP10 Rabbit pAb
Gene ID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-139 of coronavirus NSP10 (YP_009725306.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTD1
Protein Weight	794kDa
Shipping	Ice bag