

Serpina10 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uk55

Item Number: CM0026739



Introduction

SERPINA10 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for coagulation research on protein Z-dependent protease inhibitor. 51kDa.

[Learn More](#)

Parameter	Value
Product Name	SERPINA10 Rabbit pAb
Remarks/Alias	PZI; ZPI; SERPINA10
Species	Human
Gene ID	51156
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UK55
Protein Weight	51kDa
Shipping	Ice bag