

Spire2 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Wwl2

Item Number: CM0019315



Introduction

SPIRE2 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, rat. Targets actin nucleation factor SPIRE2, involved in actin dynamics, vesicle transport, and DNA damage repair. Ideal for cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	SPIRE2 Rabbit pAb
Remarks/Alias	Spir-2; SPIRE2
Species	Human
GeneID (Human)	84501
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human SPIRE2 (NP_115827.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WWL2
Protein Weight	80kDa
Shipping	Ice bag