

Tab1 Rabbit Mab For Wb, If/Icc, Elisa - Q15750

Item Number: CM0004309



Introduction

High-quality TAB1 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat TAB1. Ideal for innate immunity and MAPK signaling research.

[Learn More](#)

Parameter	Value
Product Name	TAB1 Rabbit mAb
Remarks/Alias	3'-Tab1; MAP3K7IP1; TAB1
Species	Human
Gene ID	10454
Immunogen	Synthetic peptide within amino acids 405-504 of human TAB1 (Q15750).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15750
Protein Weight	55kDa
Shipping	Ice bag