


Trimethyl-Histone H3 (K9) Rabbit Polyclonal Antibody - Q16695/P68431 | Wb, Chip, Ihc, Elisa

Item Number: CM0000395



Introduction

High-quality rabbit polyclonal antibody targeting trimethylation at lysine 9 of histone H3 (Q16695/P68431). Suitable for Western blot, ChIP, immunohistochemistry, and ELISA applications in multiple species.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K9 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; TriMethyl-Histone H3-K9
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide - A synthetic trimethylated peptide around K9 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag