


Trmt112 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ui30

Item Number: CM0014873



Introduction

Rabbit polyclonal antibody against human TRMT112, a multifunctional methyltransferase subunit involved in rRNA, tRNA, and protein methylation. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for epigenetic and translation research.

[Learn More](#)

Parameter	Value
Product Name	TRMT112 Rabbit pAb
Remarks/Alias	TRMT112; HSPC152; HSPC170; hTrm112; TRMT11-2; TRMT112
Species	Human
GeneID (human)	51504
GeneID	51504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human TRMT112 (NP_057488.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UI30
Protein Weight	14kDa
Shipping	Ice bag