


Tsfm Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P43897

Item Number: CM0003934



Introduction

TSFM Rabbit Monoclonal Antibody for detecting mitochondrial elongation factor Ts in WB, IF/ICC, IP, ELISA. Cross-reacts with human, mouse, rat. Suitable for mitochondrial translation research.

[Learn More](#)

Parameter	Value
Product Name	TSFM Rabbit mAb
Remarks/Alias	EFTS; EFTSMT; TSFM
Species	Human
GeneID (Human)	10102
GeneID	10102
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TSFM (P43897).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43897
Protein Weight	35kDa
Shipping	Ice bag