

Usp18 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Umw8

Item Number: CM0022495

Introduction

High-quality USP18 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat USP18. Ideal for studying ISG15-specific protease function in interferon signaling and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	USP18 Rabbit pAb
Remarks/Alias	ISG43; UBP43; PTORCH2; USP18
Species	Human
Gene ID (Human)	11274
Gene ID	11274
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human USP18 (NP_059110.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UMW8
Protein Weight	43kDa
Shipping	Ice bag