

Usp7/Hausp Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q93009

Item Number: CM0009663



Introduction

Rabbit monoclonal antibody targeting human USP7/HAUSP (Q93009). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for p53 pathway, DNA repair, and deubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	USP7/HAUSP Rabbit mAb
Remarks/Alias	TEF1; HAUSP; HAFOUS; DEL16P13.2; C16DELp13.2; USP7/HAUSP
Species	Human
GeneID (Human)	7874
GeneID	7874
Immunogen	Recombinant protein (recombinant fusion protein containing a sequence corresponding to amino acids 100–200 of human USP7/HAUSP (Q93009))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93009
Protein Weight	128kDa
Shipping	Ice bag