

Vamp8 Rabbit Mab For Wb, Ihc-P, Elisa - Q9Bv40

Item Number: CM0012122



Introduction

Anti-VAMP8 rabbit monoclonal antibody validated for WB, IHC-P, ELISA; reacts with human, mouse, and rat. Suitable for studying SNARE-mediated membrane fusion, autophagy, and exocytosis.

[Learn More](#)

Parameter	Value
Product Name	VAMP8 Rabbit mAb
Remarks/Alias	EDB; VAMP-8; VAMP8
Species	Human
GeneID (Human)	8673
GeneID	8673
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VAMP8 (NP_003752.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV40
Protein Weight	11kDa
Shipping	Ice bag