

Wnt3A Rabbit Pab For Wb, Elisa - P56704

Item Number: CM0015993



Introduction

WNT3A Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting human and mouse Wnt-3a. Ideal for canonical Wnt signaling, developmental biology, and stem cell self-renewal research.

[Learn More](#)

Parameter	Value
Product Name	WNT3A Rabbit pAb
Remarks/Alias	WNT3A
Species	Human
GeneID	89780
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human WNT3A (NP_149122.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P56704
Protein Weight	39kDa
Shipping	Ice bag