


Yes1 Rabbit Pab For Wb, If/Icc, Elisa - P07947

Item Number: CM0014380



Introduction

YES1 Rabbit Polyclonal Antibody (pAb) validated for WB, IF/ICC, ELISA. Recognizes Human and Mouse Tyrosine-protein kinase Yes (c-yes), a regulator of cell growth, adhesion, and migration. Ideal for cancer research and signal transduction studies.

[Learn More](#)

Parameter	Value
Product Name	YES1 Rabbit pAb
Remarks/Alias	Yes; c-yes; HsT441; P61-YES; YES1
Species	Human
GeneID (human)	7525
GeneID	7525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human YES1 (NP_005424.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07947
Protein Weight	61kDa
Shipping	Ice bag