

Znf701 Rabbit Pab For Wb, Elisa - Q9Nv72

Item Number: CM0019138

Introduction

ZNF701 Rabbit Polyclonal Antibody for Western blot (WB) and ELISA applications. Detects human ZNF701, a zinc finger protein involved in transcriptional regulation. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	ZNF701 Rabbit pAb
Remarks/Alias	ZNF701
Species	Human
GeneID (Human)	55762
GeneID	55762
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human ZNF701
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NV72
Protein Weight	61kDa
Shipping	Ice bag