

2'-O-Methylcytidine/Cm Rabbit Polyclonal Antibody For Dot Blot And Elisa Applications

Item Number: CM0000020

Introduction

Rabbit polyclonal antibody targeting 2'-O-methylcytidine (Cm), a modified nucleoside. Suitable for dot blot (DB) and ELISA applications with species-independent cross reactivity. Ideal for epitranscriptomics research.

[Learn More](#)

Parameter	Value
Product Name	2'-O-methylcytidine/ Cm Rabbit pAb
Remarks/Alias	Cm; 2'-O-methylcytidine; 2'-O-methylcytidine/ Cm
Species	Human
Immunogen	Chemical compounds
Source	Rabbit
Category	Other Modified Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag