

2'-O-Methylguanosine (Gm) Rabbit Polyclonal Antibody For Dot Blot And Elisa

Item Number: CM0000024



Introduction

Rabbit polyclonal antibody targeting 2'-O-methylguanosine (Gm) modification. Suitable for dot blot (DB) and ELISA applications. Species-independent reactivity, ideal for epitranscriptomics research.

[Learn More](#)

Parameter	Value
Product Name	2'-O-Methylguanosine(Gm) Rabbit pAb
Species	Human
Immunogen	Chemical compounds Gm
Source	Rabbit
Category	Other Modified Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag