

2'-O-Methyluridine (Um) Rabbit Polyclonal Antibody For Elisa And Dot Blot - Um

Item Number: CM0023851

Introduction

Rabbit polyclonal antibody targeting 2'-O-methyluridine (Um), a prevalent RNA modification. Suitable for ELISA and dot blot (DB) with species-independent reactivity, supporting epitranscriptomic research and diagnostic assay development.

[Learn More](#)

Parameter	Detail
Product Name	2'-O-methyluridine(Um) Rabbit pAb
Remarks/Alias	Um; 2'-O-methyluridine; 2'-O-methyluridine(Um)
Species	Human
Immunogen	Chemical compounds Um
Source	Rabbit
Category	Polyclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag