

4-Hydroxynonenal (4-Hne) Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, And Flow Cytometry (Intracellular)

Item Number: CM0000068

Introduction

The 4-Hydroxynonenal (4-HNE) Rabbit Monoclonal Antibody is a species-independent tool for detecting lipid peroxidation adducts. Validated for WB, IF/ICC, IF-P, and intracellular flow cytometry, ideal for oxidative stress research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	4-Hydroxynonenal Rabbit mAb
Remarks/Alias	4-HNE
Immunogen	Chemical compounds
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, IF/ICC, IF-P, FC (intra)
Cross Reactivity	Species independent
Shipping	Ice bag