

5-Methylcytidine/M5C Rabbit Polyclonal Antibody - M5C

Item Number: CM0020217

Introduction

Species-independent rabbit polyclonal antibody against 5-methylcytidine (m5C) for dot blot (DB) and ELISA applications. Suitable for DNA/RNA methylation studies in diagnostic and research settings.

[Learn More](#)

Parameter	Specification
Product Name	5-methylcytidine/m5C Rabbit pAb
Immunogen	Chemical compounds m5C
Source	Rabbit
Category	Polyclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag