

5-Methylcytosine (5Mc) Rabbit Pab - 5Mc

Item Number: CM0021209

Introduction

High-quality rabbit polyclonal antibody targeting 5-Methylcytosine (5mC), applicable in dot blot (DB) and ELISA. Species-independent cross reactivity enables broad epigenetic research applications.

[Learn More](#)

Parameter	Details
Product Name	5-Methylcytosine (5mC) Rabbit pAb
Immunogen	Chemical compounds 5mC
Source	Rabbit
Category	Polyclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag