

5T4 Rabbit Monoclonal Antibody For Wb And Elisa - Q13641

Item Number: CM0003947



Introduction

High-specificity rabbit monoclonal antibody against human 5T4 (TPBG), validated for Western blot and ELISA. Ideal for cancer biology, Wnt signaling, and trophoblast research - Q13641.

[Learn More](#)

Parameter	Value
Product Name	5T4 Rabbit mAb
Remarks/Alias	5T4; M6P1; 5T4AG; WAIF1
Species	Human
Gene ID (Human)	7162
Gene ID	7162
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 321-420 of human 5T4 (Q13641).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13641
Protein Weight	46kDa
Shipping	Ice bag