

A-Actin-1 (Acta1) Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P68133

Item Number: CM0005580



Introduction

High-quality rabbit monoclonal antibody targeting human α -Actin-1 (ACTA1), validated for WB, IF-P, and ELISA. Reacts with human, mouse, rat. Ideal for skeletal muscle research and myopathy studies.

[Learn More](#)

Parameter	Value
Product Name	α -Actin-1 (ACTA1) Rabbit mAb
Remarks/Alias	ACTA; ASMA; CFTD; MPFD; NEM1; NEM2; NEM3; SHPM; CFTD1; CFTDM; CMYP2A; CMYP2B; CMYP2C; α -Actin-1 (ACTA1)
Species	Human
GeneID (human)	58
GeneID	58
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human α -Actin-1 (ACTA1) (P68133).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68133
Protein Weight	42kDa
Shipping	Ice bag