


A-Synuclein Rabbit Pab - P37840

Item Number: CM0025876



Introduction

Anti- α -Synuclein rabbit polyclonal antibody validated for WB, IF/ICC, IF-P, ELISA. Detects human, mouse, rat α -Synuclein, a neuronal protein regulating synaptic vesicle trafficking and neurotransmitter release.

[Learn More](#)

Parameter	Specification
Product Name	α -Synuclein Rabbit pAb
Remarks/Alias	PD1; NACP; PARK1; PARK4; α -Synuclein
Species	Human
GeneID (Human)	6622
GeneID	6622
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-140 of human alpha-Synuclein (NP_000336.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37840
Protein Weight	14kDa
Shipping	Ice bag