


Abcc8 Rabbit Pab For Wb, Elisa - Q09428

Item Number: CM0026455



Introduction

ABCC8 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, recognizing human ABCC8/SUR1 with cross-reactivity to mouse. Essential for insulin secretion and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	ABCC8 Rabbit pAb
Remarks/Alias	HI; SUR; HHF1; MRP8; PHH1; SUR1; ABC36; HRINS; PNDM3; TNDM2; SUR1delta2; ABCC8
Species	Human
GeneID (Human)	6833
GeneID	6833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-1000 of human ABCC8 (NP_000343.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q09428
Protein Weight	177kDa
Shipping	Ice bag