

Abcf3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nuq8

Item Number: CM0019305

Introduction

ABCF3 Rabbit pAb is a polyclonal antibody suitable for WB, IHC-P, and ELISA applications. It detects human, mouse, and rat ABCF3, a protein involved in antiviral defense against flaviviruses. Ideal for immunodetection research.

[Learn More](#)

Parameter	Details
Product Name	ABCF3 Rabbit pAb
Remarks/Alias	EST201864; ABCF3
Species	Human
GeneID (Human)	55324
GeneID	55324
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ABCF3 (NP_060828.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NUQ8
Protein Weight	80kDa
Shipping	Ice bag