

Abhd12 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N2K0

Item Number: CM0023224

Introduction

ABHD12 Rabbit Polyclonal Antibody validated for WB, ELISA. Targets human ABHD12, cross-reacts with rat. Suitable for studying lysophosphatidylserine lipase function in neurological and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ABHD12 Rabbit pAb
Remarks/Alias	PHARC; ABHD12A; BEM46L2; hABHD12; C20orf22; dj965G21.2; ABHD12
Species	Human
GeneID (Human)	26090
GeneID	26090
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-250 of human ABHD12 (NP_056415.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8N2K0
Protein Weight	45kDa
Shipping	Ice bag