

Abhd5 Rabbit Pab - Q8Wts1

Item Number: CM0028867

Introduction

High-quality rabbit polyclonal antibody against ABHD5 (CGI-58, Q8WTS1), validated for WB and ELISA, cross-reacting with human and mouse. Supports lipid metabolism, triacylglycerol regulation, and keratinocyte differentiation research.

[Learn More](#)

Parameter	Value
Product Name	ABHD5 Rabbit pAb
Remarks/Alias	CGI58; IECN2; NCIE2; ABHD5
Species	Human
GeneID (human)	51099
GeneID	51099
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human ABHD5 (NP_057090.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WTS1
Protein Weight	39kDa
Shipping	Ice bag