

Acbd5 Rabbit Pab - Q5T8D3

Item Number: CM0012864



Introduction

High-quality rabbit polyclonal antibody against human ACBD5 (Q5T8D3). Validated for Western blot and ELISA; suitable for pexophagy, lipid metabolism, and peroxisome research. Cross-reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	ACBD5 Rabbit pAb
Remarks/Alias	RDLKD; ACBD5
Species	Human
GeneID (Human)	91452
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-440 of human ACBD5 (NP_663736.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5T8D3
Protein Weight	60kDa
Shipping	Ice bag