

Acetyl-Alpha Tubulin K40 Rabbit Monoclonal Antibody (Q71U36) For Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000175



Introduction

Recombinant rabbit monoclonal antibody targeting acetylated K40 of human alpha-tubulin (TUBA1A, Q71U36). Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for cytoskeleton and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-alpha Tubulin -K40 Rabbit mAb
Remarks/Alias	LIS3; TUBA3; B-ALPHA-1; Acetyl-alpha Tubulin -K40
Species	Human
GeneID (Human)	7846
GeneID	7846
Immunogen	Synthetic peptide A synthetic acetylated peptide around K40 of human alpha Tubulin(NP_001257328.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q71U36
Protein Weight	50 kDa
Shipping	Ice bag