

Acetyl-Alpha Tubulin (K40) Rabbit Pab (Q71U36) For Wb, If-P, Elisa

Item Number: CM0000108



Introduction

Rabbit polyclonal antibody specific for acetylated lysine 40 of human alpha-tubulin (TUBA1A, Q71U36). Validated for WB, IF-P, and ELISA; cross-reacts with mouse and rat tubulin. Ideal for acetylation and microtubule dynamics studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-alpha Tubulin -K40 Rabbit pAb
Remarks/Alias	LIS3; TUBA3; B-ALPHA-1; Acetyl-alpha Tubulin -K40
Species	Human
GeneID	7846
Immunogen	Synthetic acetylated peptide around K40 of human alpha Tubulin
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q71U36
Protein Weight	50kDa
Shipping	Ice bag