


Acetyl-Brd4 (K332) Rabbit Polyclonal Antibody (Swiss O60885) | Wb, Elisa - Brd4 Acetylation Detection

Item Number: CM0000116



Introduction

Acetyl-BRD4-K332 Rabbit pAb (O60885) specifically detects acetylated BRD4 at lysine 332 in human and mouse samples. Validated for Western Blot and ELISA, this antibody is ideal for studying epigenetic regulation and transcription mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-BRD4-K332 Rabbit pAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; Acetyl-BRD4-K332
Species	Human
GeneID (Human)	23476
GeneID	23476
Immunogen	Synthetic Peptide A synthetic acetylated peptide around K332 of human BRD4 (NP_490597.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60885
Protein Weight	152kDa
Shipping	Ice bag