

Acetyl-Dnmt1 (K1127/K1129/K1131/K1133) Rabbit Polyclonal Antibody (P26358) For Ihc-P, If/Icc, Elisa

Item Number: CM0000115



Introduction

Rabbit polyclonal antibody targeting acetylated DNMT1 at K1127, K1129, K1131, K1133 (UniProt P26358). Validated for IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. 183 kDa. Ideal for DNA methylation and epigenetics studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-DNMT1-K1127/K1129/K1131/K1133 Rabbit pAb
Remarks/Alias	AIM; DNMT; MCMT; CXXC9; HSN1E; ADCADN; m.Hsal; Acetyl-DNMT1-K1127/K1129/K1131/K1133
Species	Human
GeneID (human)	1786
GeneID	1786
Immunogen	A synthetic acetylated peptide around K1127 & K1129 & K1131 & K1133 of human DNMT1 (NP_001124295.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26358
Protein Weight	183kDa
Shipping	Ice bag