


Acetyl-Histone H1.4-K26 Rabbit Polyclonal Antibody (P10412) | IHC-P, IF/ICC, Elisa

Item Number: CM0000152



Introduction

Acetyl-Histone H1.4-K26 Rabbit pAb (P10412) validated for IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H1.4-K26 Rabbit pAb
Remarks/Alias	H1E; H1.4; H1F4; RMNS; H1s-4; HIST1H1E; dj221C16.5; Acetyl-Histone H1.4-K26
Species	Human
GeneID (human)	3008
GeneID	3008
Immunogen	Synthetic peptide: A synthetic acetylated peptide around K26 of human Histone H14 (NP_005312.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10412
Protein Weight	22kDa
Shipping	Ice bag