


Acetyl-Histone H2A (K5) Rabbit Polyclonal Antibody (P04908/P0C0S8) For Wb, If/Icc, Elisa

Item Number: CM0000167



Introduction

Rabbit polyclonal antibody targeting acetylated lysine 5 (K5ac) of human Histone H2A (SWISS P04908/P0C0S8). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat, and predicted broad species. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H2A-K5 Rabbit pAb
Remarks/Alias	H2A/c; H2AFC; H2AC11; H2AC15; H2AC16; H2AC17; HIST1H2AI; Acetyl-Histone H2A-K5
Species	Human
GeneID (human)	H2A
GeneID	3012/8329
Immunogen	Synthetic Peptide: A synthetic acetylated peptide around K5 of human Histone H2A (NP_003508.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P04908/P0C0S8
Protein Weight	14kDa
Shipping	Ice bag