


Acetyl-Histone H2B-K12 Rabbit Pab (P62807/Q16778) For Wb, Ihc-P, Chip

Item Number: CM0000110



Introduction

Rabbit polyclonal antibody specific for acetylated histone H2B at lysine 12 (P62807/Q16778). Validated for WB, IHC-P, IF/ICC, ELISA, ChIP. Reacts with human, mouse, rat and other species. Ideal for chromatin and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H2B-K12 Rabbit pAb
Remarks/Alias	H2B; H2BE; H2BQ; GL105; H2B.1; H2BFQ; H2BGL105; H2B-GL105; HIST2H2BE; Acetyl-Histone H2B-K12
Species	Human
GeneID (human)	3017/8349
GeneID	H2B
Immunogen	A synthetic acetylated peptide around K12 of human Histone H2B (NP_003519.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62807/Q16778
Protein Weight	14kDa
Shipping	Ice bag