

Acetyl-Histone H2B-K15 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa | P62807/Q16778

Item Number: CM0000100



Introduction

Rabbit polyclonal antibody targeting acetylated lysine 15 of Histone H2B (SWISS: P62807/Q16778). Validated for Western blot, IHC-P, IF/ICC, and ELISA. Predicted to react with human, mouse, rat, and other species.

[Learn More](#)

Parameter	Description
Product Name	Acetyl-Histone H2B-K15 Rabbit pAb
Species	Human
GeneID (Human)	H2B
GeneID	3017/8349
Immunogen	Synthetic Peptide: A synthetic acetylated peptide around K15 of human Histone H2B (NP_003519.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62807/Q16778
Shipping	Ice bag