


Acetyl-Histone H2B-K20 Rabbit Polyclonal Antibody (P62807/Q16778) For Wb, Elisa, Chip

Item Number: CM0000162



Introduction

Acetyl-Histone H2B-K20 Rabbit pAb specifically detects acetylated lysine 20 of human Histone H2B (SWISS P62807/Q16778). Validated for WB, ELISA, and ChIP in human, mouse, and rat samples. Ideal for epigenetic research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H2B-K20 Rabbit pAb
Remarks/Alias	H2B; H2BE; H2BQ; GL105; H2B.1; H2BFQ; H2BGL105; H2B-GL105; HIST2H2BE; Acetyl-Histone H2B-K20
Species	Human
GeneID (Human)	H2B
GeneID	3017/8349
Immunogen	Synthetic Peptide: A synthetic acetylated peptide around K20 of human Histone H2B (NP_003519.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P62807/Q16778
Protein Weight	14kDa
Shipping	Ice bag