


## Acetyl-Histone H2B-K5 Rabbit Polyclonal Antibody (P62807/Q16778) For Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000127



### Introduction

Acetyl-Histone H2B-K5 Rabbit pAb from CamelBio, specific for acetylated lysine 5 of human histone H2B (SWISS P62807/Q16778), suitable for Western blot, IHC, IF/ICC, and ELISA applications. Predicted wide cross-reactivity with human, mouse, rat, and other species for epigenetic research.

[Learn More](#)

| Parameter        | Value                                                                                          |
|------------------|------------------------------------------------------------------------------------------------|
| Product Name     | Acetyl-Histone H2B-K5 Rabbit pAb                                                               |
| Species          | Human                                                                                          |
| GeneID (human)   | H2B                                                                                            |
| GeneID           | 3017/8349                                                                                      |
| Immunogen        | Synthetic Peptide: A synthetic acetylated peptide around K5 of human Histone H2B (NP_003519.1) |
| Source           | Rabbit                                                                                         |
| Category         | Acetylated Antibodies                                                                          |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                       |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                                                |
| SWISS            | P62807/Q16778                                                                                  |
| Shipping         | Ice bag                                                                                        |