

Acetyl-Histone H3-K14 Rabbit Monoclonal Antibody (Q16695 / P68431) For Wb, Db, Chip, Elisa

Item Number: CM0000096



Introduction

Acetyl-Histone H3-K14 Rabbit mAb validated for WB, DB, ChIP, ELISA. Recognizes acetylated K14 of human Histone H3 (Q16695/P68431) with cross-reactivity to mouse and rat. Suitable for chromatin and epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K14 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	A synthetic acetylated peptide around K14 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag