

Acetyl-Histone H3-K14 Rabbit Monoclonal Antibody (Q16695/P68431) For Wb, Chip, Elisa

Item Number: CM0000182



Introduction

High-quality rabbit monoclonal antibody against acetylated Histone H3 at lysine 14, validated for WB, IF-P, IHC-P, ChIP, and ELISA. Cross-reacts with human, mouse, and rat. Essential for epigenetic research and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K14 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K14
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic acetylated peptide around K14 of human H3 (NP_003520.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF-P, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag