


Acetyl-Histone H3 (K14) Rabbit Polyclonal Antibody | Q16695 / P68431 | Wb, If/Icc, Elisa

Item Number: CM0000155



Introduction

Rabbit polyclonal antibody targeting acetylated Histone H3 at lysine 14 (K14). Validated for Western blot, immunofluorescence, and ELISA. Detects human, mouse, rat, and other species (wide range predicted). Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K14 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K14
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic acetylated peptide around K14 of human H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag