


Acetyl-Histone H3-K18 (P68431/Q16695) Rabbit Mab | Western Blot, Chip, Elisa

Item Number: CM0000130



Introduction

High-specificity rabbit monoclonal antibody targeting acetylated K18 of Histone H3 (P68431/Q16695). Validated for WB, DB, IHC-P, ELISA, ChIP, ChIP-seq, CUT&Tag with wide species reactivity. Ideal for epigenetic research and chromatin profiling.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K18 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; H3C1; Acetyl-Histone H3-K18
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic acetylated peptide around K18 of human Histone H3 (P68431).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, ELISA, ChIP, ChIP-seq, CUT&Tag
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag