

Acetyl-Histone H3-K18 (Q16695/P68431) Rabbit Polyclonal Antibody For Wb, Ihc-P, Chip And Elisa

Item Number: CM0000180



Introduction

Highly specific rabbit polyclonal antibody targeting acetylated lysine 18 of human Histone H3 (SWISS: Q16695/P68431). Validated for WB, DB, IHC-P, IF/ICC, ELISA, and ChIP. Ideal for epigenetic research, chromatin studies, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K18 Rabbit pAb
Remarks/Alias	H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FJ; H3C10; H3C11; HIST1H3J; Acetyl-Histone H3-K18
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic acetylated peptide around K18 of human H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag