

Acetyl-Histone H3-K18 Rabbit Monoclonal Antibody (Q16695/P68431) | Validated For Wb, Ihc, Chip, Elisa

Item Number: CM0000187



Introduction

Rabbit monoclonal antibody recognizing acetylated lysine 18 (K18) of human Histone H3 (Swiss-Prot Q16695/P68431). Validated for WB, IHC-P, IF/ICC, ChIP, and ELISA. Predicted wide cross-reactivity with human, mouse, rat, and other species.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K18 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Acetyl-Histone H3-K18
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	A synthetic acetylated peptide around K18 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	16kDa
Shipping	Ice bag